

**Table 4**

Male and Female LV Wall thickness (mm)

Age	Sex	Exposure Duration	Control	BPA 2.5	BPA 25	BPA 250	BPA 2500	BPA 25000	EE2 0.05	EE2 0.5
PND 21	Male	–	1.85 ± 0.17 (10)	1.81 ± 0.15 (8)	1.74 ± 0.11 (8)	1.90 ± 0.32 (9)	1.89 ± 0.14 (10)	1.99 ± 0.32 (10)	1.88 ± 0.22 (10)	1.95 ± 0.32 (10)
	Female	–	1.83 ± 0.18 (10)	1.91 ± 0.21 (10)	1.85 ± 0.10 (10)	1.84 ± 0.16 (10)	1.89 ± 0.12 (10)	1.87 ± 0.20 (10)	1.74 ± 0.25 (10)	1.72 ± 0.16 (10)
PND 90	Male	Stop	2.57 ± 0.44 (10)	2.64 ± 0.46 (10)	2.46 ± 0.22 (8)	2.49 ± 0.23 (10)	2.44 ± 0.31 (9)	2.48 ± 0.28 (10)	2.57 ± 0.37 (10)	2.37 ± 0.27 (10)
	Male	Continuous	2.58 ± 0.23 (10)	2.37 ± 0.28 (9)	2.38 ± 0.28 (8)	2.38 ± 0.20 (10)	2.52 ± 0.45 (8)	2.45 ± 0.22 (10)	2.49 ± 0.27 (10)	2.43 ± 0.26 (10)
	Female	Stop	2.60 ± 0.29 (10)	<b>2.30 ± 0.22 (8)<sup>a</sup></b>	2.49 ± 0.20 (10)	2.47 ± 0.28 (10)	2.49 ± 0.20 (8)	2.49 ± 0.22 (10)	2.55 ± 0.23 (10)	2.59 ± 0.30 (10)
	Female	Continuous	2.36 ± 0.22 (10)	2.53 ± 0.39 (10)	2.28 ± 0.26 (10)	2.40 ± 0.17 (10)	2.46 ± 0.30 (9)	2.22 ± 0.22 (10)	2.47 ± 0.26 (10)	2.51 ± 0.18 (10)
6 Month	Male	Stop	2.81 ± 0.43 (10)	2.68 ± 0.20 (10)	2.66 ± 0.38 (9)	2.72 ± 0.28 (10)	2.53 ± 0.29 (10)	2.75 ± 0.25 (9)	2.46 ± 0.18 (10)	2.43 ± 0.22 (10)
	Male	Continuous	2.54 ± 0.30 (10)	2.73 ± 0.43 (9)	2.62 ± 0.37 (9)	<b>2.95 ± 0.47 (10)<sup>c</sup></b>	2.69 ± 0.38 (8)	2.65 ± 0.32 (10)	2.61 ± 0.25 (10)	2.60 ± 0.36 (9)
	Female	Stop	2.44 ± 0.19 (10)	2.34 ± 0.18 (10)	2.40 ± 0.25 (10)	2.68 ± 0.27 (10)	2.73 ± 0.41 (10)	2.49 ± 0.30 (10)	2.44 ± 0.33 (10)	2.62 ± 0.24 (10)
	Female	Continuous	2.28 ± 0.36 (10)	2.39 ± 0.27 (10)	2.45 ± 0.33 (10)	2.38 ± 0.26 (10)	2.51 ± 0.24 (10)	2.49 ± 0.23 (10)	<b>2.63 ± 0.35 (10)<sup>b</sup></b>	2.44 ± 0.27 (9)

Values represent group mean ± SD; the number of samples analyzed is shown in parenthesis.

<sup>a</sup> p = .049, d = 1.19.<sup>b</sup> p = .048 d = 1.03.<sup>c</sup> p = .08 d = 1.10